

# **METALLURGICAL MICROSCOPES**



## **KERN Pictograms**



360° rotatable microscope head



Monocular Microscope For the inspection with one eve



Binocular Microscope For the inspection with both eyes



Trinocular Microscope

For the inspection with both eyes and the additional option for the connection of a camera



Abbe Condenser

With high numerical aperture for the concentration and the focusing of light



Halogen illumination

For pictures bright and rich in contrast



**LED illumination** 

Cold, energy saving and especially long-life illumination



Incident illumination

For non-transparent objects



Transmitting illumination

For transparent objects



Fluorescence illumination

For stereomicroscopes



Fluorescence illumination

for compound microscopes With 100 W mercury lamp and filter



Fluorescence illumination for compound microscopes With 3W LED illumination and filter



Phase contrast unit

For a higher contrast



Darkfield condenser/unit

For a higher contrast due to indirect illumination



Polarising unit

To polarise the light



POLAR

Infinity system

Infinity corrected optical system



Zoom magnification

For stereomicroscopes



Parallel optical system

For stereomicroscopes, enables fatigue-proof working



Integrated scale In the eveniece



SD card

For data storage



USB 2.0 digital camera

For direct transmitting of the picture to a PC



USB 3.0 digital camera

For direct transmitting of the picture to a PC



WLAN data interface:

For transmitting of the picture to a mobile display device



**HDMI** digital camera

For direct transmitting of the picture to a display device



PC software

To transfer the measurements from the device to a PC.



Automatic temperature compesation

For measurements between 10 °C and 30 °C



Protection against dust and water splashes IPxx

The type of protection is shown by the pictogram.



**Battery operation** 

Ready for battery operation. The battery type is specified for each device.



Battery operation rechargable

Prepared for a rechargable battery operation



Mains adapter

230V/50Hz in standard version for EU. On request GB, AUS or USA version.



Power supply

Integrated in microscope. 230V/50Hz standard EU. More standards e.g. GB, AUS or USA on request.



Package shipment

The time required to manufacture the product internally is shown in days in the pictogram.



C-Mount Adapter for the connection of a

camera to a trinocular microscope

Frames per second

High (Super) Wide Field (Eyepiece with H(S)WF

high eye point for wearers of glasses)

LWD Long Working Distance **Numerical Aperture** N.A.

SLR camera Single-Lens Reflex camera

SWF Super Wide Field (Field number at

least Ø 23 mm for 10× eyepiece)

W.D. Working Distance

WF Wide Field (Field number up to

Ø 22 mm for 10× eyepiece)

## Why you should choose a KERN microscope now!

For 175 years, KERN & SOHN has been synonymous with high precision weighing and measuring technology. This claim is the driving force for the development of our microscope and refractometer ranges.

Thanks to consistent customer focus paired with smart ideas and the latest available technology we are proud to be suppliers of high-quality, durable top microscopes and refractometers, which help you to be as efficient as possible in your daily work.

When developing our microscopes we have concentrated on the very best optical quality and have used only high-quality glass and the latest technologies to achieve this. The high-quality Philips halogen and modern LED illumination produce razor-sharp images with high contrast and which will impress you with their brilliant true-colour display - you must have noticed this yourself.

#### Your advantages:

- all mechanical parts have been designed for a long service life
- special attention has been given to the ergonomy of our microscopes, as this allows the user to work for several hours in a comfortable position which does not cause fatigue
- · our microscopes are fully-equipped and can be used immediately
- Highlight for 2020: die KERN camera software you will be amazed at how user-friendly and intuitive it is, a high-quality tablet camera as well as a comprehensive range of calibration services for refractometers
- · and much more...

Use our practical "Check list for microscopes and refractometers", which may help you to quickly determine specifications for the future instrument. Together with our KERN product specialist you can choose the right product for you.

If there is no suitable product in the standard range, for example, then we will of course configure an individual microscope for you.

Our aim is to develop a market-driven product solution, so with our microscope and refractometer range, the saying holds true: good quality at a competitive price! This is what we stand for and work towards, every day!

With our current 2020 product range you can benefit from improved quality and a clear reduction in price, which we have been able to achieve through more efficient production methods and increased global sales of our microscopes and refractometers and of course we pass this straight on to you

Do you have any questions about our range of microscopes and

Your KERN customer consultants are available at any time to help you further.

I hope that you enjoy working efficiently with our KERN Optics products.

Albert Sauter, Managing Director

## Your advantages

#### fast

- · 24 hour dispatch service for products in stock - order today, on its way tomorrow
- · Sales & service hotline from 8:00 am to 5:00 pm

#### reliable

- Up to 3 years warranty
- · Certified QM system **DIN EN ISO 9001**

#### versatile

- · One-stop shopping: from microscope through to refractometer - everything from one supplier
- Quick as a flash, find the product you want with the "Quick-Finder" at www.kern-sohn.com



## Order hotline

+49 7433 9933-0



# Service hotline

+49 7433 9933-199



## Online shop

www.kern-sohn.com



#### E-mail order

info@kern-sohn.com



#### **Fax Order**

+49 7433 9933-146

will assist you



## Our team of consultants

from Monday to Friday from 8:00 am to 5:00 pm



#### www.kern-sohn.com

Information on current product availability, product data sheets, user instructions, useful knowledge, technical glossary and much for you to download, practical topic areas, which will guide you to the right product in your industry as well as a microscopes and refractometers search engine.

# 2 Metallurgical microscopes







Illumination unit with filter disc



Stage and objectives

#### LAB LINE MET

## The metallurgical reflected light microscope for material testing and surface testing, as well as quality assurance in industry

#### **Features**

- The KERN OKM is an excellent metallurgical reflected light microscope, e.g. for surface quality testing of raw materials and finished products in industry
- The strong, continuously dimmable 30 W halogen reflected illumination unit (Philips) ensures excellent, high-contrast images
- · The illumination unit with an integrated 5-slot filter wheel for blue, green, yellow, grey and blank means that you can quickly change the colour filter for different contrast views
- · A large mechanical stage for reflected illumination applications is configured as standard. The coarse and fine focusing knob on both sides guarantees optimal adjustment and focusing of your sample

- A simple polarising unit (analyser and polariser) is included with delivery
- · A large selection of different eyepieces, objectives and a polarising unit are also
- A protective dust cover, eye cups, as well as multi-lingual user instructions are included in the scope of delivery
- A C-mount adapter is required to connect a camera to the trinocular version. You can select this adapter from the following model outfit list
- · Please find detailed information in the following model outfit list

#### Scope of application

· Metallurgy, material testing, quality assurance

#### Applications/Samples

· Opaque and thick samples, workpieces (surfaces, fold lines, coatings)

### **Technical data**

- · Infinity optical system
- · Quadplex nosepiece
- Siedentopf 30° inclined/360° rotatable
- · Diopter adjustment: One-sided
- · Overall dimensions W×D×H 440×200×460 mm
- · Net weight basic configuration approx. 8 kg



Model	Standard configuration					
KERN	Tube	Eyepiece	Objective quality	Objectives	Illumination	ex works €
OKM 172*	Binocular	HWF 10×/Ø 18 mm	Infinity Plan	5×/10×/ LWD 20×/	30 W Halogen (incident)	1790,-
OKM 173	Trinocular	HWF 10×/Ø 18 mm	Infinity Plan	LWD40×	30 W Halogen (incident)	1890,-

Model outfit		Mode	I KERN	Order number	Price/piece excl. of VAT
		OKM 172	OKM 173		ex works €
	HWF 10×/Ø 18 mm	✓	✓	OBB-A1403	45,-
	HWF 10×/ø 18 mm (reticule 0,1 mm) (non-adjustable)	✓	✓	OBB-A1349	60,-
Eyepieces (23,2 mm)	WF 5×/Ø 20 mm	0	0	OBB-A1355	55,-
(==,=,	WF 12,5×/ø 14 mm	0	0	OBB-A1353	75,-
	WF 16×/ø 13 mm	0	0	OBB-A1354	45,-
	5×/0,11 W.D. 12,10 mm	✓	✓	OBB-A1268	105,-
Infinity Plan achromatic	10×/0,25 W.D. 4,75 mm	✓	✓	OBB-A1244	195,-
objectives (no cover glass)	20×/0,40 (spring-loaded) W.D. 2,14 mm	0	0	OBB-A1251	265,-
(110 dover glads)	40×/0,65 (spring-loaded) W.D. 0,45 mm	0	0	OBB-A1258	290,-
Infinity	20×/0,40 W.D. 8,35 mm	✓	✓	OBB-A1252	290,-
Plan achromatic objectives	40×/0,65 W.D. 3,90 mm	✓	✓	OBB-A1259	405,-
(no cover glass) for long working	50×/0,70 (spring-loaded) W.D. 1,95 mm	0	0	OBB-A1266	450,-
distance	80×/0,80 (spring-loaded) W.D. 0,85 mm	0	0	OBB-A1271	500,-
Binocular tube	Siedentopf 30° inclined/360° rotatable     Interpupillary distance 50 – 75 mm     Diopter adjustment: One-sided	✓	0	OBB-A1130	
Trinocular tube	Siedentopf 30° inclined/360° rotatable     Interpupillary distance 50 – 75 mm     Light distribution 80:20     Diopter adjustment: One-sided	0	~	OBB-A1346	
Mechanical stage	Stage size W×D 200×140 mm     Travel 76×52 mm     Coaxial coarse and fine focusing knobs	<b>✓</b>	<b>*</b>		
Illumination	30 W Halogen spare bulb (incident)	✓	✓	OBB-A1372	35,-
Reflected illumination unit	5-filter unit (Blue, Green, Yellow, Grey, Empty)	✓	✓		
	Polarising unit (Incl. analyser and polariser slide)	✓	<b>✓</b>		
0.14	1×		0	OBB-A1514	100,-
C-Mount	0,5× (focus adjustable)		0	OBB-A1515	160,-

✓ = Included with delivery

O = Option





Stage OKN



Stage OKO



Illumination unit

#### **PROFESSIONAL LINE MET**

## The fully-equipped reflected and transmitted light microscope for numerous applications in metallurgy

#### **Features**

- · The KERN OKN and OKO series are professional, versatile, metallurgical microscopes, which are used in testing metals and analysing surfaces
- · You can choose between two reflected illumination units: A 50 W halogen reflected illumination unit or a premium illumination unit with stunning 100 W reflected illumination for powerful performance
- A height-adjustable 1,25 Abbe condenser which can be centred as well as a field diaphragm are available for the transmitted light variants (KERN OKO series), for complete professional Koehler illumination
- The KERN OKO transmitted illumination variant is fitted with an open, mechanical stage, as standard. On the other hand, the KERN OKN reflected illumination variant has a closed, mechanical stage, as standard

- · A simple polarising unit (analyser and polariser) is included with delivery
- · A large selection of accessories, such as, for example, a Butterfly tube, eyepieces and further objectives are available for longer working distances
- A protective dust cover, eye cups, as well as multi-lingual user instructions are included in the scope of delivery
- A C-Mount adapter is required to connect a camera. You can select this adapter from the following model outfit list
- · Please find detailed information in the following model outfit list

#### Scope of application

· Metallurgy, material testing, quality assurance

#### Applications/Samples

· Opaque and thick samples, workpieces (surfaces, fold lines, coatings)

### **Technical data**

- · Infinity optical system
- Quintuple nosepiece
- Siedentopf 30° inclined/360° rotatable
- · Diopter adjustment: Both-sided
- · Overall dimensions W×D×H 550×200×460 mm
- · Net weight basic configuration approx. 14,5 kg

#### STANDARD























#### Model Standard configuration **Price** excl. of VAT Tube ex works Eyepiece Objective quality Objectives Illumination KERN 5×/10×/LWD20×/ OKN 175\* Trinocular WF 10×/Ø 18 mm Infinity Plan 50 W Halogen (incident) 2230,-LWD40× 5×/10×/LWD20× 50 W Halogen (incident) + **OKO 176\*** Trinocular WF 10×/Ø 18 mm Infinity Plan 2770,-LWD40×/100× 20 W (transmitted) 5×/10×/LWD20×/ **OKN 177** 100 W Halogen (incident) Trinocular WF 10×/Ø 18 mm Infinity Plan 2450,-LWD40> 5×/10×/LWD20×/ 100 W Halogen (incident) + **OKO 178** Trinocular 2990,-WF 10×/Ø 18 mm Infinity Plan LWD40×/100× 20 W (transmitted)

\*ONLY WHILE STOCKS LAST

Model outfit			Mode	KERN	Order number	Price/piece excl. of VAT	
			OKO OKN 176 177		OKO 178		ex ci. of val ex works €
	WF 10×/Ø 18 mm	175 ✓ <b>√</b>	11	44	11	OBB-A1347	40,-
Eyepieces (23,2 mm)	WF 10×/Ø 18 mm (reticule 0,1 mm) (adjustable)		1	✓	✓	OBB-A1350	90,-
	WF 5×/Ø 20 mm	00	00	00	00	OBB-A1355	55,-
(20,2 11111)	WF12,5×/Ø 14 mm	00	00	00	00	OBB-A1353	75,-
	WF 16×/Ø 13 mm	00	00	00	00	OBB-A1354	45,-
	5×/0,11 W.D. 6,73 mm	✓	✓	✓	✓	OBB-A1268	105,-
	10×/0,25 W.D. 4,19 mm	✓	✓	✓	✓	OBB-A1244	195,-
Infinity Plan objectives	20×/0,40 (spring-loaded) W.D. 2,14 mm	0	0	0	0	OBB-A1251	265,-
(no cover glass)	40×/0,65 (spring-loaded) W.D. 0,45 mm	0	0	0	0	OBB-A1258	290,-
	100×/1,25 (oil) (spring-loaded) W.D. 0,12 mm		✓		✓	OBB-A1241	315,-
Infinity	20×/0,40 W.D. 8,35 mm	✓	✓	✓	✓	OBB-A1252	290,-
Plan objectives	40×/0,65 W.D. 3,90 mm	✓	✓	✓	✓	OBB-A1259	405,-
(no cover glass) for long working	50×/0,70 (spring-loaded) W.D. 1,95 mm	0	0	0	0	OBB-A1266	450,-
distance	80×/0,80 (spring-loaded) W.D. 0,85 mm	0	0	0	0	OBB-A1271	500,-
Trinocular tube	Siedentopf 30° inclined/360° rotatable     Interpupillary distance 50 – 75mm     Light distribution 100:0     Diopter adjustment: Both-sided	<b>✓</b>	~	~	~		
	Butterfly 30° inclined/360° rotatable     Interpupillary distance 50 – 75 mm     Light distribution 100:0     Diopter adjustment: Both-sided	0	0	0	0	OBB-A1382	520,-
Mechanical stage for reflection	Stage size W×D 200×140 mm Travel 78×55 mm Stage fast lowering unit Stage Up-Down moving range: max. 50 mm	~		~			
Mechanical stage for transmitted Ilumination	Stage size W×D 175×145 mm     Travel 78×55 mm     Coaxial coarse and fine focusing knobs		✓		<b>✓</b>		
Stage plate	Plate for sample placement	✓	1	✓	✓		
Glass plate	Glass plate		0		0	OBB-A1378	50,-
Reflected illumination unit	Polarising unit (Incl. analyser, polariser and blue filter slide)	~	✓	✓	✓		
Condenser	Abbe N.A. 1,25 (aperture diaphragm)		1		✓	OBB-A1380	100,-
	20 W Halogen spare bulb (transmitted)		1		✓	OBB-A1370	25,-
llumination	50 W Halogen spare bulb (incident)	✓	1			OBB-A1207	19,-
	100 W Halogen spare bulb (incident)			✓	✓	OBB-A1377	35,-
Polarising unit	for transmitted illumination		✓		✓	OBB-A1470	135,-
	Blue		1		✓		
Colour filters	Green		0		0	OBB-A1188	19,-
or transmitted Ilumination	Yellow		0		0	OBB-A1165	19,-
	Grey		0		0	OBB-A1183	19,-
	1×	0	0	0	0	OBB-A1140	100,-
C-Mount	0,57× (focus adjustable)	0	0	0	0	OBB-A1136	160,-

✓ = Included with delivery

O = Option





Specimen stage and illumination unit



Analyser/Polariser

#### **LAB LINE MET**

## The inverted metallurgical microscope for professional applications

#### **Features**

- The KERN OLM range is part of the range of inverted microscopes and stands out through its design which is ergonomic, robust and extremely stable. This range, with its large working distance is, for example, particularly suitable for surface quality testing of raw materials and finished products in industry
- · Strong and continuously adjustable 50W halogen illumination unit ensures the optimum illumination of the materials to be tested
- · As standard, the OLM range is fitted with a trinocular eyepiece tube
- · A simple polarising unit (analyser and polariser) is included with delivery

- · A large mechanical stage is included with delivery as standard. The coarse and fine focusing knob on both sides guarantees optimal adjustment and focusing
- Further options such as, for example, a large selection of objectives can be integrated as accessories
- · A dust cover as well as user instructions are included with the delivery
- · Please find detailed information in the following model outfit list

#### Scope of application

· Metallurgy, material testing, quality assurance

#### Applications/Samples

• Opaque and thick samples, workpieces (surfaces, fold lines, coatings)

### **Technical data**

- · Infinity optical system
- Quintuple nosepiece
- Siedentopf 30° inclined
- · Diopter adjustment: Both-sided
- · Overall dimensions W×D×H 271×379×747 mm
- · Net weight approx. 12,5 kg

















ŧ.	
٧	1 DAY

Model	Standard configuration					
KERN	Tube	Eyepiece	Objective quality	Objectives	Illumination	ex works €
OLM 171	Trinocular	HWF 10×/Ø 22 mm	Infinity Plan	LWD5×/LWD10×/ LWD20×/LWD50×	50 W Halogen (incident)	3190,-

Model outfit		Model KERN	Order number	Price/piece excl. of VAT
		OLM 171		ex works €
Evepieces	HWF 10×/ø 22 mm (adjustable)	✓	OBB-A1491	90,-
(30 mm)	HWF 10×/ø 22 mm (reticule 0,1 mm) (adjustable)	✓	OBB-A1523	140,-
	5×/0,13 W.D. 24,23 mm	✓	OBB-A1525	125,-
Infinity	10×/0,25 W.D. 18,48 mm	✓	OBB-A1526	185,-
Plan achromatic	20×/0,40 W.D. 8,35 mm	✓	OBB-A1527	245,-
<b>objectives</b> for long working	50×/0,70 (spring-loaded) W.D. 1,95 mm	✓	OBB-A1528	320,-
distance	80×/0,80 (spring-loaded) W.D. 0,85 mm	0	OBB-A1530	480,-
	100×/0,90 (dry) W.D. 1,0 mm	0	OBB-A1531	720,-
Trinocular tube	30° inclined     Interpupillary distance 48-76 mm     Light distribution 100:0     Diopter adjustment: Both-sided	<b>√</b>		
Mechanical stage	Stage size W×D 210×180 mm     Travel 50×50 mm     Coaxial coarse and fine focusing knobs	<b>✓</b>		
Illumination	50 W Halogen spare bulb (incident)	✓	OBB-A1207	19,-
Reflected illumination unit	Polarising unit (Incl. analyser, polariser and colour filter slide)	✓		
	Blue	✓	OBB-A1510	19,-
Colour filters for transmitted illumination	Green	0	OBB-A1511	19,-
	Yellow	0	OBB-A1512	19,-
	Grey	0	OBB-A1513	19,-
0.14	0,5×	0	OBB-A1515	185,-
C-Mount	1×	0	OBB-A1514	120,-

<sup>✓ =</sup> Included with delivery

Printed in Germany by KERN & SOHN GmbH z-co-gb-kp-20201

+49 7433 9933-0 Tel. Fax +49 7433 9933-146

info@kern-sohn.com www.kern-sohn.com







inclination scale out of iron, a foundation stone for weighing for the balance



1844



1880



1923



1980

ousts mechanical devices



1994





2002

certification in accordance with DIN EN ISO 9001:2000



2007

of medical products (DIN EN 13485 and



2008

Authorisation for initial verification by the



2009

Approval for the manufacture and sale of height rods (DIN EN 13485



2012

Verification point for non-automatic balances



Expansion of the product range to include optical instruments (microscopes



Inauguration of Ziegelei 2.0 with computer-controlled



Come with KERN into the



#### 2019

Significant anniversary year 250 years of balance manufacture in the Sauter family-owned company



2020

Construction of Ziegelei 3.0,